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1. Carson, J.L. et al. 1987. Rapid, easy and economical screening tests for antibodies to human immunodeficiency virus. Lancet ii:361-362.
2. Van de Perre, P.D. et al. 1988. Comparison of six serological assays for human immunodeficiency virus antibody detection in developing countries. J. Clin. Microbiol. 26:552-556.

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MedMira Laboratories Inc.
加拿大麦美华公司
155 Chain Lake Drive, Suite 1
Halifax, Nova Scotia, Canada
B3S 1B3
Telephone: 902 450 1588
Fax: 902 450 1580
www.medmira.com (English/英文)
www.medmira.cn (Chinese/中文)

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Rapid HIV Test Kit

This package insert must be read carefully prior to using the product. Accuracy of test results cannot be guaranteed if there are deviations from the enclosed instructions.

PACKAGING FORMAT

- Single test
- Box of 30 tests

INTENDED USE

Rapid HIV test kit can detect antibodies to HIV-1 and/or HIV-2 in whole blood. It can be used in blood donation, general public health, emergency, blood transfusion and remote areas without proper equipment.

BIOLOGICAL PRINCIPLES OF THE TEST

Human immunodeficiency virus (HIV) causes acquired immunodeficiency syndrome (AIDS). Infection with HIV-1 and/or HIV-2 results in the production of the corresponding HIV antibodies. Therefore, tests that detect HIV-1/HIV-2 antibodies provide the means to identify HIV infected individuals.^{1,2} Rapid HIV Test Kit utilizes synthetic peptides corresponding to conserved regions of the HIV-1 and HIV-2 *env* genomes for the qualitative detection of antibodies to HIV-1 group M & O and/or HIV-2. Sample is added to the reaction zone. A conjugate cap contains the indicator reagent, producing a colorimetric reaction that can be easily visualized and interpreted. An integrated vertical red procedural control line appears to indicate that the reagent and test components are functioning properly and that the test procedure has been performed correctly. When the conjugate binds to the HIV antibodies, a horizontal test line also appears. In the absence of HIV antibodies, only the procedural control line appears.

CONTENTS

Single test

- 1 Mylar Pouch
- 1 test cartridge with attached Blue colour Specimen Filtration Unit (For Whole Blood Test)
- 1 Instant Gold™ Cap (Green colour)
- 1 disposable pipette
- 1 Universal Buffer dropper bottle (White colour Cap)
- 1 vial of Diluent Buffer Solution (Pink colour Cap)

Box of 30 tests

- 30 Mylar Pouches containing:
 - 1 test cartridge with attached Blue colour Specimen Filtration Unit (For Whole Blood Test)
 - 1 Instant Gold™ Cap (Green colour)
 - 1 disposable pipette
 - 1 Universal Buffer dropper bottle (White colour Cap)
 - 1 vial of Diluent Buffer Solution (Pink colour Cap)

Note: Each test provides one set of lancet and alcohol swab

STORAGE CONDITION AND EXPIRY DATE

The Rapid HIV Test Kit should be stored in a cool dry area protected from direct sunlight, heat and moisture. Do not freeze, store at room temperature 2-30 °C
Expiry Date: 24 months

EQUIPMENT REQUIRED

No additional equipments are required.

SPECIMEN COLLECTION

For Whole Blood

- Using the alcohol swab provided, clean your Index finger. Allow the finger to dry thoroughly.
- Pull the white plastic tip off the lancet and place it against your clean site of your Index finger. Puncture the skin by pressing the handle all the way. Hold the finger downward and apply gentle pressure beside the point of puncture.
- Squeeze the outside of the disposable pipette and place it in the center of the blood drop. Slowly release channel of the pipette. Apply **one (1)** drop of blood into the uncapped Diluent Buffer vial (Pink colour cap). Gently mix by tapping the bottom of the vial. Immediately begin the **TESTING PROCEDURE**.

INSTRUCTIONS

Please read the following items carefully before performing the test.

- Read this instruction sheet completely and carefully prior to use of the Rapid HIV Test Kit. If the directions are not followed exactly, inaccurate test results may occur.
- Do not use after the expiration date printed on the outside box or on the Mylar pouch
- Store in a dry place at 2-30°C
- Keep out of reach of children
- For *in vitro* diagnostic use only. Not to be taken internally
- Do not open the Mylar pouch until you are ready to start the test
- Dispose with care for all test kit materials after use

DIRECTIONS FOR USE

GENERAL PREPARATION

- Remove the Mylar pouch from the box.
- Using the notched corners, tear open the Mylar pouch. Place all components on a clean, flat surface.
- Ensure the blue Specimen Filtration Unit is placed firmly in the well of the white plastic Test Cartridge.

TESTING PROCEDURE

- All solutions must be completely absorbed into the test membrane before proceeding to the next step in the testing procedure.
- Once the assay has been started, all subsequent steps should be completed without interruption.
- Perform the test on a flat work surface.
- Read the test results immediately. Failure to do so may result in inaccurate test results.

FOR WHOLE BLOOD SPECIMEN

- Apply six (6) drops from the MedMira Universal Buffer vial (White colour cap) to the center of the Blue Specimen Filtration Unit. Allow the buffer to absorb completely. Pour the entire contents of the Diluent Buffer vial (Pink Colour cap) into the center of the blue Specimen Filtration Unit. Ensure that the solution has completely absorbed into the blue Specimen Filtration Unit before proceeding to the next step.
- Remove the blue Specimen Filtration Unit by gently twisting the handle in a circular and upwards motion. Place the instant Gold™ Cap (Green colour) loosely into the well of the test cartridge. Apply **twelve (12)** drops of MedMira Universal Buffer to the center of the Instant Gold™ Cap.
- Allow the buffer to absorb completely. Remove the Instant Gold™ Cap and **Read Result**.

QUALITY CONTROL

To insure assay validity, a procedural control is incorporate in the device. If there is no procedural control line under the C, and even if the line appears adjacent to the T, the test is invalid and should be repeated with a new Rapid HIV Test Kit.

HOW TO READ THE TEST RESULTS

You Probably Are Not Exposed to HIV-1 and/or HIV-2

The presence of one red line in the control region (C) means you probably are not exposed to HIV-1 and/or HIV-2. If the symptoms still persist or there is reason to be concern, repeat the procedure with another Rapid HIV Test Kit or consult your doctor.

You might be exposed to HIV-1 or HIV-2

The presence of two red lines in both the test region (T) and control region (C) means you might be exposed to HIV-1 and/or HIV-2. One line may be lighter than the other, and they do not have to match. It means that HIV-1 and/or HIV-2 antibodies are probably present in your blood. Visit your doctor as soon as possible.

Invalid Result

The result is invalid if no red line appears in the control region (C), even if a line appears in the test region (T). Also, the presence of a broken line under the (C) indicates that there has been a problem, either with the test device or the specimen, during the Testing Procedure. You should repeat the procedure with a new Rapid HIV Test Kit. If the problem persists, contact your local distributor.

LIMITATIONS OF THE PROCEDURE

- The procedure must be followed as directed to obtain accurate results.
- The intensity of the test line does not necessarily correlate to the titre of antibody in the specimen.
- No serological test can provide absolute assurance that a specimen does not contain the low level of HIV- 1/2 antibodies present at an early stage of infection.
- A negative result therefore does not exclude the possibility of exposure to or infection with HIV-1/2 viruses.
- All positive specimens must be confirmed by supplemental testing.

PERFORMANCE CHARACTERISTICS

Sensitivity

The sensitivity of the Rapid HIV Test Kit for the detection of antibodies to HIV-1 group M subtypes was determined to be 100%. The Rapid HIV Test Kit showed 100% sensitivity when 10 confirmed HIV-1 group O specimens were tested. The sensitivity of the Rapid HIV Test Kit for HIV-2 antibodies was determined to be 100% when 99 Western blot confirmed HIV-2 positive specimens were tested (Table 1). The overall sensitivity of the Rapid HIV Test Kit was 99.93% when 2036 Western blot confirmed HIV-1 positive specimens (different clades) were tested (Table 2).

Table 1 Sensitivity of Rapid HIV Test Kit with HIV-1 group M, group O and HIV-2

Specimen (group, clade)	Number of Samples	Rapid HIV Test Kit Sensitivity (%)
HIV-1 (M,A)	8	100
HIV-1 (M,B)	58	100
HIV-1 (M,C)	19	100
HIV-1 (M,D)	4	100
HIV-1 (M,E)	3	100
HIV-1 (M,F)	3	100
HIV-1 (M,G)	1	100
HIV-1 (O)	10	100
HIV-2	99	100

Table 2 Sensitivity of Rapid HIV Test Kit with prevalent HIV-1 clades by region

Country	Prevalence	Number of Samples	Rapid HIV Test Kit Sensitivity (%)
Canada	B	836	99.76
India	C,A,B	20	100
Peru	B,F	20	100
Tanzania	C,A,D	14	100
Thailand	E,B	20	100
Trinidad	B	20	100
United States	B	1105	99.73

Specificity

The overall specificity of the Rapid HIV Test Kit was 99.31% when 15,420 HIV-1/HIV-2 negative specimens were tested.

Reproducibility

There were neither intra- or inter-observer variations nor intra- or inter-lot variations, and the reproducibility of the Rapid HIV Test Kit was found to be 100%.

Comparability with a FDA Approved Rapid HIV Test

The Rapid HIV Test Kit showed 100% concordance with the FDA approved SUDS rapid HIV test (Table 3).

Table 3 Comparative Study between Rapid HIV Test Kit and SUDS

SUDS result	Rapid HIV Test Kit Result		Total
	+	-	
+	50	0	50
-	0	150	150
total	50	150	200

Interference Study

Interference study was performed on 400 specimens and the results indicated that EDTA, heparin , sodium citrate and elevated analytes AST (range 96 – 1764 U/L), ALT (136 – 2296 U/L), ALP (297 – 550 U/L), total bilirubin (34 – 65 µmol/L), LD (350 – 600 U/L) and Uric acid (600 – 940 µmol/L) did not interfere with the test.

PRODUCT WARRANTY

MedMira Laboratories Inc. guarantees the quality of this product if stored and used as stipulated. Any component of the kit found to be defective shall be replaced free of charge upon return of the defective product. MedMira Laboratories Inc. disclaims any implied warranty of merchantability or fitness for a particular purpose, and in no event shall the MedMira Laboratories Inc. be liable for consequent damage.

LITERATURE CITED

- Carson, J.L. et al. 1987. Rapid, easy and economical screening tests for antibodies to human immunodeficiency virus. *Lancet* ii : 361-362.
- Van de Perre, PD. et al. 1988. Comparison of six serological assays for human immunodeficiency virus antibody detection in developing countries. *J. Clin. Microbiol.* 26: 552- 556.

MedMira Laboratories Inc. Distributor
Suite 1, 155 Chain Lake Drive
Halifax, Nova Scotia. B3S 1B3
Canada
Telephone: +1-902-450-1588
Fax: +1-902-450-1580
Email: sales@medmira.com
www.medmira.com
April 2003 CINS0140 Rev 1/2

CFDA Registration Number: 20143401076
Registration Specification: YZB/CAN 0594-2014
Product insert approval date and revision date: February 24, 2014