

multipl

rapid HBc/HIV/HCV antibody test

Rapid co-infection detection: three diseases, one test.



Multipl Rapid HBc/HIV/HCV Antibody Test.

Multipl Rapid HBc/HIV/HCV Antibody Test simultaneously detects antibodies for HIV-1/2, hepatitis B and hepatitis C in a single test.

It is the only rapid test for these diseases that uses MedMira's patented Rapid Vertical Flow Technology™. In three simple steps, you get an easy-to-read answer, going from sample collection to results in minutes.

Testing for these three diseases at the same time can lead to speedier, more effective, treatment. Patients do not have to wait for answers, consultations can happen in a timely manner and costs are reduced.

Features

- Multiple diseases tested on a single cartridge
- Simple test procedure
- Instant, easy-to-interpret results
- Built-in procedural and reagent control line
- Works with whole blood, serum or plasma
- Room-temperature storage
- No timers or specialized equipment needed
- Product-specific training not required

Applications

These three diseases share common risk factors and similar transmission routes. Oftentimes, patients and medical professionals are not aware that co-infection exists. Yet co-infection is common, reaching rates of 30–35% for HIV and HCV. By testing for these diseases at the same time, patient outcomes can be improved since more complete therapeutic action plans can be initiated.

Here are just a few of the applications where Mutiplo HBc/HIV/HCV is improving human health.

Public Health Programs

Voluntary Testing & Counseling Centers

Pharmacies & Physician Offices

Immigration Screening

Military Medical Units


Convenience Care Clinics

Occupational Exposures


Blood, Organ & Tissue Banks



One test. More answers.



Non-Reactive Result



Reactive Result

Ordering Information

POC for fingerstick whole blood

Cat. No.: 815311004551

This format is ideally suited for:

- Point-of-care patient testing
- Settings requiring a portable and all-in-one test
- Mobile clinics and field-testing

LAB+ for venipuncture whole blood/serum/plasma

Cat. No.: 815311004575

This format is ideally suited for:

- Testing serum, plasma, and/or whole blood specimens
- Maximum flexibility in specimen type and testing environment

LAB S/P for serum/plasma

Cat. No.: 8153211004582

This format is ideally suited for:

- Laboratories and hospitals
- Testing serum or plasma specimens
- Batch testing

For more information on Multiplo HbC/HIV/HCV, please contact MedMira's Sales & Marketing team at sales@medmira.com or +1 902 450 1588

Technology & Performance

Multiplo HbC/HIV/HCV takes advantage of MedMira's patented Rapid Vertical Flow Technology. This technology has delivered excellent clinical performance in rigorous evaluations and inspections.

	HBc	HIV	HCV
Sensitivity	96.2%	99.8%	99.1%
Specificity	100%	99.7%	99.7%

About MedMira

MedMira's patented Rapid Vertical Flow Technology delivers fast, accurate results in a range of applications. Our products help care providers and their patients know more, instantly. MedMira rapid tests are manufactured in a cGMP facility with a quality management system registered to ISO 9001:2008 and ISO 13485:2003. Learn more at medmira.com